

	1 x 1 min 70% ethanol
- distilled water	2 x 1 min
- washing buffer (e.g.PBS/BSA)	2 x 2min
- blocking solution	background reactions
- primary antibody	up to 24 hours depending on antibody concentration
- washing buffer	3 x 5 min
- GOD labeled secondary antibody (other detection principle, e.g. biotin-avidin)	30 min
- PBS	3 x 5 min
- enzyme substrate, chromogen	20-60 min (variable incubation)
- PBS	2 x 2 min
- distilled water	2 x 2 min
- counterstain	optional
- distilled water	2 x 2 min
- mounting	glycerol/gelatine or other mounting

Two-step enzyme technique (GOD as label and HRP as indicator)

Enzyme substrate (GOD)	Cytochemical feature
<p>0.05% Diaminobenzidine (DAB), 1% D-glucose containing 0.1-1 mg peroxidase (HRP) per mL solution</p> <ul style="list-style-type: none"> • DAB stock solution: 50 mg DAB are dissolved in 10 mL distilled water, aliquoted and stored at -20°C (stable for years) • β-D-Glucose stock solution: 1000 mg D-glucose dissolved in 10 mL distilled water (for a final 1% glucose concentration), aliquoted and stored at 4°C • Peroxidase (HRP RZ 3) stock solution: 10 mg HRP dissolved in 10 mL 0.1 M phosphate buffer pH 7.3 (for 0.1 mg peroxidase per mL working solution) • Working solution: 1 mL DAB stock <i>plus</i> 3 mL 0.1 M phosphate buffer pH 6.8 <i>plus</i> 1 mL D-glucose stock solution <i>plus</i> 1 mL peroxidase stock solution <i>plus</i> 4 mL 0.1 M phosphate buffer pH 6.8 	<p>Brown color</p> <p>Enzyme substrate is useful in combination with other substrates for double stainings</p>
<p>General procedure for immunohistological staining (paraffin sections):</p> <ul style="list-style-type: none"> - Deparaffinize in xylene 2 x 5 min - rehydrate with ethanol 2 x 5 min absolute ethanol 1 x 1 min 95% ethanol 	

	1 x 1 min 70% ethanol
- distilled water	2 x 1 min
- washing buffer (e.g. PBS/BSA)	2 x 2min
- blocking solution	background reactions
- primary antibody	up to 24 hours depending on antibody concentration
- washing buffer	3 x 5 min
- GOD labeled secondary antibody (other detection principle, e.g. biotin-avidin)	30 min
- washing buffer	3 x 5 min
- enzyme substrate, chromogen	10-30 min (variable incubation)
- washing buffer, distilled water	2 x 2 min, each
- postfixation with 0.1% OsO ₄	2-3 min (optional)
- counterstain	haematoxylin or other (optional)
- dehydration, mounting	ascending series of ethanol, xylene, mounting medium

Enzyme substrate (GOD)	Cytochemical feature
<p>0.05% Diaminobenzidine (DAB), 1% D-glucose, 0.025% nickel ammonium sulfate containing 0.1-1 mg peroxidase (HRP) per mL solution</p> <ul style="list-style-type: none"> • DAB stock solution: 50 mg DAB are dissolved in 10 mL distilled water, aliquoted and stored at -20°C (stable for years) • β-D-Glucose stock solution: 1000 mg D-glucose dissolved in 10 mL distilled water (for a final 1% glucose concentration), aliquoted and stored at 4°C • Peroxidase (HRP RZ 3) solution: 10 mg HRP dissolved in 10 mL 0.1 M phosphate buffer pH 7.3 (for 0.1 mg peroxidase per mL working solution) • Nickel ammonium sulfate stock solution: 0.1 g Ni(NH₄)₂ (SO₄)₂ are dissolved in 10 mL distilled water • Working solution: 1 mL DAB stock <i>plus</i> 3 mL 0.1 M phosphate buffer pH 6.8 <i>plus</i> 1 mL D-glucose stock solution <i>plus</i> 1 mL peroxidase stock solution <i>plus</i> 3.75 mL 0.1 M phosphate buffer (0.01 M pH 6.8) <i>plus</i> 250 µL nickel ammonium sulfate stock 	<p>Black color</p> <p>Enzyme substrate is used for enhanced GOD stainings and useful in combination with other substrates for double stainings</p>
<p>General procedure for immunohistological staining (paraffin sections):</p>	
<p>- Deparaffinize in xylene</p> <p>- rehydrate with ethanol</p>	<p>2 x 5 min</p> <p>2 x 5 min absolute ethanol</p> <p>1 x 1 min 95% ethanol</p>

	1 x 1 min 70% ethanol
– distilled water	2 x 1 min
– washing buffer (e.g.PBS/BSA)	2 x 2min
– blocking solution	background reactions
– primary antibody	up to 24 hours depending on antibody concentration
– washing buffer	3 x 5 min
– GOD labeled secondary antibody (other detection principle, e.g. biotin-avidin)	30 min
– washing buffer	3 x 5 min
– enzyme substrate, chromogen	10-30 min at room temperature or 1-2 hours at 4°C (optimal incubation time and temperature depend on experimental conditions)
– washing buffer, distilled water	2 x 2 min, each
– counterstain	haematoxylin or other (optional)
– dehydration, mounting	ascending series of ethanol, xylene, mounting medium

Two-step enzyme technique (GOD as label and HRP as indicator) in the coimmobilized GOD-HRP method

In this enzyme-amplification technique, H₂O₂ generated by the reaction of GOD with D-glucose serves as substrate for HRP in the oxidation of chromogens such as 3,3'-diaminobenzidine (DAB). Coimmobilization of both *marker enzyme* (GOD) and *secondary system amplifier* (HRP) in close proximity proved to be of great advantage in immunohistology because the unfavorable properties of soluble HRP in substrate mixtures are avoided. This method appears to be better than the use of soluble HRP in the substrate mixture (KUHLMANN WD and PESCHKE P, 1986).

The coupled GOD-HRP principle was developed for the use of GOD as marker enzyme in immunohistology because no such endogenous enzyme activity is known in mammalian cells and to profit from the general usefulness of DAB as cytochemical stain. Hence, the coupled GOD-HRP principle can be superior to conventional immunoperoxidase labeling for the localization of biomolecules in tissues which are rich in endogenous peroxidase activities.

Enzyme substrate (GOD) (coimmobilized GOD-HRP method)	Cytochemical feature
<p>0.05% Diaminobenzidine (DAB), 0.2-1% β -D-glucose 0.04% ammonium chloride (optional)</p> <ul style="list-style-type: none"> • 0.5% Cobalt acetate for pretreatment is optional: 0.5 g cobalt acetate dissolved in 100 mL Tris buffer pH 7.6 • Ammonium chloride stock solution is optional: 0.4 g NH₄Cl dissolved in 100 mL distilled water • DAB stock solution: 50 mg DAB are dissolved in 10 mL distilled water, 	<p>Brown color, dark blue/black color</p> <p>GOD-HRP method as alternative to immunoperoxidase staining, useful in combination with other markers and substrates</p>

References for further readings

Nachlas MM *et al.* (1957)
Reiss J (1966)
Kuhlmann WD (1970)
Kuhlmann WD and Avrameas S (1971)
Altman FP (1976)
Clark CA *et al.* (1982)
Gay H *et al.* (1984)
Kuhlmann WD (1984)
Kuhlmann WD and Peschke P (1986)

Full citation of publications is given in chapter *References*

link: <https://www.kuhlmann-biomed.de/wp-content/uploads/2020/12/References.pdf>