

| | |
|---|--|
| – blocking solution | endogenous peroxidases |
| – blocking solution | background reactions |
| – primary antibody | up to 24 hours depending on antibody concentration |
| – washing buffer | 3 x 5 min |
| – HRP labeled secondary antibody (other detection principle, e.g. biotin-avidin) | 30 min |
| – PBS | 3 x 5 min |
| – 0.2 M Tris-HCl buffer pH 7.2-7.4 | rinse |
| – enzyme substrate, chromogen | 10-20 min (variable incubation) |
| – PBS | 2 x 2 min |
| – distilled water | 2 x 2 min, |
| – postfixation with 0.1% OsO ₄ | 1-3 min (optional) |
| – counterstain | haematoxylin or other (optional) |
| – dehydration, mounting | ascending series of ethanol, xylene, mounting medium |

Note: intensity of HRP staining can be enhanced at an acid pH or at a neutral pH by the addition of imidazole to the DAB/H₂O₂ substrate mixture. The nitrogenous compound was found to raise the enzyme activity of HRP in test tube assays by about 30%.

Furthermore, the addition of cobalt chloride and nickel ammonium sulfate will intensify the DAB reaction product

Aminoethylcarbazole technique

| Enzyme substrate (HRP) | Cytochemical feature |
|--|---|
| <p>0.05% Aminoethylcarbazole* (AEC) and 0.01% H₂O₂</p> <ul style="list-style-type: none"> • AEC stock solution: 0.1 g AEC are dissolved in 10 mL DMF** aliquoted and stored at 4°C (stable for years) • H₂O₂ stock solution: 100 µL 30% H₂O₂ dissolved in 3.0 mL distilled water • Working solution: 250 µL AEC stock dissolved in 4.7 mL acetate buffer (0.05 M pH 5.5) <i>plus</i> 50.0 µL H₂O₂ stock <p>* 3-Amino-9-ethylcarbazole ** N,N-Dimethylformamide</p> | <p>Red color</p> <p>The substrate is also useful in combination with other substrates for double stainings</p> |
| <p>General procedure for immunohistological staining (paraffin sections):</p> <ul style="list-style-type: none"> – Deparaffinize in xylene 2 x 5 min – rehydrate with ethanol 2 x 5 min absolute ethanol 1 x 1 min 95% ethanol 1 x 1 min 70% ethanol | |

| | |
|---|---|
| - distilled water | 2 x 1 min |
| - washing buffer (e.g.PBS/BSA) | 2 x 2min |
| - blocking solution | endogenous peroxidases |
| - blocking solution | background reactions |
| - primary antibody | up to 24 hours depending on antibody concentration |
| - washing buffer | 3 x 5 min |
| - HRP labeled secondary antibody (other detection principle, e.g. biotin-avidin) | 30 min |
| - PBS | 3 x 5 min |
| - 0.05 M acetate buffer pH 5.5 | rinse |
| - enzyme substrate, chromogen | 5-15 min at 37°C in the dark staining intensity is controlled under the microscope |
| - distilled water | 3 x 2 min |
| - counterstain | haematoxylin or other (optional) |
| - mounting | glycerol/gelatine or other mounting medium |

Chloro-naphthol technique

| Enzyme substrate (HRP) | Cytochemical feature |
|---|--|
| <p>0.03% Chloro-naphthol (CN)* and 0.01% H₂O₂</p> <ul style="list-style-type: none"> • CN stock solution: 0.3 g CN are dissolved in 10 mL absolute ethanol aliquoted and stored at -20°C (stable for months) • H₂O₂ stock solution: 100 µL 30% H₂O₂ dissolved in 3.0 mL distilled water • Working solution: 100 µL CN stock dissolved in 9.8 mL Tris-HCl buffer (0.05 M pH 7.4-7.6) <i>plus</i> 100 µL H₂O₂ stock <p>Substrate solution is stable for about 60 min and must be filtered prior to use</p> <p>* 4-Chloro-1-naphthol</p> | <p>Blue color</p> <p>The substrate is also useful in combination with other substrates for double stainings</p> |
| <p>General procedure for immunohistological staining (paraffin sections):</p> <ul style="list-style-type: none"> - Deparaffinize in xylene - rehydrate with ethanol - distilled water - washing buffer (e.g.PBS/BSA) - blocking solution | |

| | |
|---|---|
| – blocking solution | background reactions |
| – primary antibody | up to 24 hours depending on antibody concentration |
| – washing buffer | 3 x 5 min |
| – HRP labeled secondary antibody (other detection principle, e.g. biotin-avidin) | 30 min |
| – PBS | 3 x 5 min |
| – 0.05 M Tris-HCl buffer pH 7.4-7.6 | rinse |
| – enzyme substrate, chromogen | 5-20 min at 37°C in the dark staining intensity is controlled under the microscope |
| – distilled water | 3 x 2 min |
| – counterstain | haematoxylin or other (optional) |
| – mounting | glycerol/gelatine or other mounting medium |

HRP-GOD coupled technique

In neuroanatomical tracing technique, modifications of the original HRP-DAB method have been described. One proposed way is the use of glucose oxidase (GOD) and D-glucose in the substrate medium by which a continuous release of hydrogen peroxide is generated instead of direct addition of hydrogen peroxide for the marker enzyme HRP (ITO H K et al., 1979; OLDFIELD B J et al., 1983; SAKANAKA M et al., 1987; SHU SY et al, 1988). It was suggested that the continuous release of hydrogen peroxide would favor the deposition of DAB reaction products around HRP molecules. This coupled enzyme method together with metal intensification was used successfully in bringing out histochemical and immunoreactive structures.

| Enzyme substrate (HRP) (HRP-GOD coupled method) | Cytochemical feature |
|--|--|
| <p>0.05% Diaminobenzidine (DAB), 0,2% β -D-glucose 0.04% ammonium chloride containing 0.3-0.5 mg glucose oxidase (GOD) per 100 mL solution</p> <ul style="list-style-type: none"> • <i>0.5% Cobalt acetate for pretreatment: 0.5 g cobalt acetate dissolved in 100 mL Tris buffer pH 7.6</i> • Ammonium chloride stock solution: 0.4 g NH₄Cl dissolved in 100 mL distilled water • DAB stock solution: 50 mg DAB are dissolved in 10 mL distilled water, aliquoted and stored at -20°C (stable for years) • β-D-Glucose stock solution: 200 mg D-glucose dissolved in 10 mL distilled water, aliquoted and stored at 4°C • Glucose oxidase (GOD, e.g. 250 U/mg) solution: 3 mg GOD dissolved in 10 mL 0.1 M phosphate | <p>Dark blue/black color</p> <p>Method to increase the sensitivity of the HRP method, useful in combination with other substrates</p> |

Graham RC and Karnovsky MJ (1966)
Nakane PK (1968)
Itoh K *et al.* (1979)
Oldfield BJ *et al.* (1983)
Sakanaka M *et al.* (1987)
Shu SY *et al.* (1988)

Full citation of publications is given in chapter *References*

link: <https://www.kuhlmann-biomed.de/wp-content/uploads/2020/12/References.pdf>

© Prof. Dr. Wolf D. Kuhlmann

10.02.2010